



# Eco-Bite

Biopsy Forceps – Single Use


## Eco-Bite

- Sharp Standard & Cups with Spike  
*Precise & clean tissue acquisition*
- Versions with Coated Metal Spiral with Positioning Bands  
*Friction minimized passage through working channels*

*Convenient evaluation of forceps position for a controlled & safe catheter advancement and removal*



**Medi-Globe®**  
*The Spirit of Care.*



## Cup Shapes



## Bronchoscopy

### Biopsy Forceps Standard

Eco-Bite – Biopsy Forceps with Metal Spiral, sterile packed				
Catalog No.	Cup Shapes	Outer Diameter	Working Length	Minimum Working Channel Diameter
GBF-31-18-120	Oval Cups	1.8 mm	120 cm	2.0 mm
GBF-32-18-120	Oval Cups with Spike	1.8 mm	120 cm	2.0 mm

### Biopsy Forceps with Coated Metal Spiral

Eco-Bite – Biopsy Forceps with Coated Metal Spiral, sterile packed				
Catalog No.	Cup Shapes	Outer Diameter	Working Length	Minimum Working Channel Diameter
GBF-41-18-120	Oval Cups	1.8 mm	120 cm	2.0 mm
GBF-42-18-120	Oval Cups with Spike	1.8 mm	120 cm	2.0 mm

## Gastroscopy & Colonoscopy

### Biopsy Forceps Standard

Eco-Bite – Biopsy Forceps with Metal Spiral, sterile packed				
Catalog No.	Cup Shapes	Outer Diameter	Working Length	Minimum Working Channel Diameter
GBF-31-23-180	Oval Cups	2.3 mm	180 cm	2.5 mm
GBF-31-23-230	Oval Cups	2.3 mm	230 cm	2.5 mm
GBF-32-23-180	Oval Cups with Spike	2.3 mm	180 cm	2.5 mm
GBF-32-23-230	Oval Cups with Spike	2.3 mm	230 cm	2.5 mm

### Biopsy Forceps with Coated Metal Spiral and Position Bands

Eco-Bite – Biopsy Forceps with Coated Metal Spiral and Position Bands, sterile packed				
Catalog No.	Cup Shapes	Outer Diameter	Working Length	Minimum Working Channel Diameter
GBF-41-23-180	Oval Cups	2.3 mm	180 cm	2.5 mm
GBF-41-23-230	Oval Cups	2.3 mm	230 cm	2.5 mm
GBF-42-23-180	Oval Cups with Spike	2.3 mm	180 cm	2.5 mm
GBF-42-23-230	Oval Cups with Spike	2.3 mm	230 cm	2.5 mm

L-GBF-SU-01-ENG Ver. 2.0 (2017-04)