

M-Flex® Blue Silicone Bougies



There’s nothing more flexible than an M-Flex®.

Our silicone bougies feature added flexibility that some users prefer. Tungsten-filled and latex-free, there’s nothing quite like the M-Flex®.

Benefits

- Ultra-flexible tapered tip (Maloney)
- Completely radiopaque from distal tip to proximal end
- American and European depth markings
- Mercury-free and latex-free
- Autoclavable for quick reprocessing

Features

- 21 sizes: 20 FR to 60 FR (in 2 FR increments)
- Choose a full set, mini set or individual pieces
- Three-year warranty on product and workmanship
- Easy disposal
- No charge trade-in program for proper disposal of expired mercury-filled bougies (with purchase of new product)




Store and transport M-Flex® bougies in our storage case which holds up to two full sets.

M-Flex® Dilatation Systems

Set	Part No.	Contents
Maloney (tapered) 21-piece full set	1206-00	21 bougies, 20–60 FR
Maloney (tapered) 10-piece mini set	1206-10	10 bougies, 36–54 FR
Hurst (blunt) 21-piece full set	1207-00	21 bougies, 20–60 FR
Hurst (blunt) 10-piece mini set	1207-10	10 bougies, 36–54 FR
Bougie Storage Case	1210-01	Heavy-duty storage for up to 46 bougies. 32"L x 7"W x 4"H




Our storage case allows easy identification of bougie size. Removing plexiglass window provides easy access.



End caps display French size.



Extra flexible tip.



Top: Maloney tip. Bottom: Hurst tip.

M-Flex® Blue Silicone Bougies

FR (French) Size	(mm)	Maloney (tapered) Tip Part No.	Hurst (blunt) Tip Part No.
20	6.6	1206-20	1207-20
22	7.4	1206-22	1207-22
24	8	1206-24	1207-24
26	8.6	1206-26	1207-26
28	9.3	1206-28	1207-28
30	10	1206-30	1207-30
32	10.7	1206-32	1207-32
34	11.2	1206-34	1207-34
36	12	1206-36	1207-36
38	12.5	1206-38	1207-38
40	13	1206-40	1207-40
42	14	1206-42	1207-42
44	14.5	1206-44	1207-44
46	15.2	1206-46	1207-46
48	16	1206-48	1207-48
50	16.5	1206-50	1207-50
52	17.2	1206-52	1207-52
54	18	1206-54	1207-54
56	18.5	1206-56	1207-56
58	19.2	1206-58	1207-58
60	20	1206-60	1207-60

Hurst Bougie

Maloney Bougie

American Depth Markings

European Depth Markings

