


Fascial Dilators Set

Fascial dilators are used for dilatation of percutaneous track over the guidewire.



- Fully radiopaque polyethylene
- Female luer lock allows heparin, contrast media and other solutions to be injected using standard luer lock syringe
- Accepts 0.038" guidewires

QTY /Box: 5

Order No.	Description	Length (cm)
FDS 6/10 FDS 8/16	Fascial Dilator Set 6, 7, 8, 9, 10 Fascial Dilator Set 8, 10, 12, 14, 16	20

Individual Fascial Dilators

QTY /Box: 5

Order No.	Description	Length (cm)
FD 6	Fascial Dilator 6F	20
FD 7	Fascial Dilator 7F	
FD 8	Fascial Dilator 8F	
FD 9	Fascial Dilator 9F	
FD 10	Fascial Dilator 10F	
FD 12	Fascial Dilator 12F	
FD 14	Fascial Dilator 14F	
FD 16	Fascial Dilator 16F	
FD 18	Fascial Dilator 18F	