

DILATATION CATHETERS

Biliary Dilatation Catheter

Used to dilate the biliary strictures, pancreatic strictures or papillary orifice.



• Length 200 cm

6

BDC

QTY/Box: 5

Order No.	Transition Length	TAPER	
		From French Size	To French Size
BDC 7/5	2.5 cm	5F	7F
BDC 10/6	2.5 cm	6F	10F

BDC -T [PTFE]

QTY/Box: 5

Order No.	Transition Length	TAPER	
		From French Size	To French Size
BDC 7/5 -T	2.5 cm	5F	7F