

Gastroscopy Biteblocks



Indications



Biteblocks used to facilitate the insertion of the instrument during endoscopy procedures.

**ITABLOC®**

Features



- Preloaded headstrap in laprene with differentiated elasticity.
- With tongue pad.
- The diameter of the central port permits the access of various dilators from ch15 to ch 60
- Available without headstrap.
- Disposable version.
- Autoclavable version
- Packaged not sterile.
- 100 pcs. per box.

References	Diameter	Material	Colour	Headstrap
MED/102	21x25 mm	polypropylene	yellow	•
MED/102/B	21x25 mm	polypropylene	yellow	•
MED/102/A	21x25 mm	polycarbonate 	transparent	•
MED/102/AB	21x25 mm	polycarbonate 	transparent	•

Specifications

- Disposable biteblock in polypropylene
- Autoclavable biteblock in polycarbonate
- Laprene headstrap with differentiated elasticity