



Oesophageal strictures close to the cardia require a specific stent design. The cardia umbrella stents have a umbrella designed distal end which significantly reduces irritations at the great curve of the stomach and prevent an in- or overgrow of tumor into the stent.

All cardia umbrella stents are preloaded on the reload delivery system. For correction of the stent position the stent can be reloaded into the delivery system during the procedure. The excellent radiopacity tantalum markers makes the stents clearly visible in the fluoroscopy.

- Proximal partially or complete silicone covering
- Radiopaque tantalum markers
- Distal end umbrella design
- Atraumatic ends
- High tensile extraction thread on proximal end
- Memory effect of nitinol material – high radial force
- Reload delivery system

Length 700 mm, OTW

Ø 8 mm / 24 French

Accepts guide wire 0.035 inch

## aixstent · cardia umbrella · partially silicone covered

REF	Length	Length Cover	Ø Middle	Ø Ends
122-24-100	100	85	24	30 / 50
122-24-120	120	105	24	30 / 50

## aixstent · cardia umbrella · complete silicone covered

REF	Length	Length Cover	Ø Middle	Ø Ends
123-24-100	100	100	24	30 / 50
123-24-120	120	120	24	30 / 50