



The colon and rectum stents are preloaded on the OTW (over-the-wire) delivery systems. Their specific end designs have the perfect fit for the anatomic and peristaltic conditions in the colon and the rectum. For correction of the stent position the stent can be reloaded into the delivery system during the procedure. The excellent radiopacity tantalum marker makes the stents clearly visible in the fluoroscopy.

- Stent without and with partially silicone covering
- Atraumatic ends
- Special tulip design prevent migration
- High tensile extraction threads on both ends
- Memory effect of nitinol material – high radial force
- Radiopaque tantalum markers
- Reload delivery system

Length 1100 mm, OTW
 Ø 8 mm / 24 French
 Accepts guide wire 0.035 inch

aixstent · colon rectum · non covered

| REF | Length | Length Cover | Ø Middle | Ø Ends |
|------------|--------|--------------|----------|--------|
| 401-30-080 | 80 | without | 30 | 36 |
| 401-30-100 | 100 | without | 30 | 36 |

aixstent · colon rectum · partially silicone covered

| REF | Length | Length Cover | Ø Middle | Ø Ends |
|------------|--------|--------------|----------|--------|
| 402-30-080 | 80 | 50 | 30 | 36 |
| 402-30-100 | 100 | 70 | 30 | 36 |