



Leufen Medical offers a huge variety of self-expandable oesophageal stents with diameters of 20 and 24 mm and lengths from 60 to 140 mm. The special end design and the strong radial force keep the stents in position and avoid stent migration.

All oesophageal stents are preloaded on the reload delivery system. For correction of the stent position the stent can be reloaded into the delivery system during the procedure. The excellent radiopacity tantalum markers makes the stents clearly visible in the fluoroscopy.

- Stent with partially or complete silicone covering
- Atraumatic ends
- High tensile extraction threads on both ends
- Memory effect of nitinol material – high radial force
- Radiopaque tantalum markers
- Reload delivery system

Length 700 mm, OTW

Ø 8 mm / 24 French

Accepts guide wire 0.035 inch

## aixstent · oesophagus · partially silicone covered

REF	Length	Length Cover	Ø Middle	Ø Ends
102-20-060	60	30	20	26
102-20-080	80	50	20	26
102-20-100	100	70	20	26
102-20-120	120	90	20	26
102-20-140	140	110	20	26
102-24-080	80	50	24	30
102-24-100	100	70	24	30
102-24-120	120	90	24	30
102-24-140	140	110	24	30

## aixstent · oesophagus · complete silicone covered

REF	Length	Length Cover	Ø Middle	Ø Ends
103-20-080	80	80	20	26
103-20-100	100	100	20	26
103-20-120	120	120	20	26
103-20-140	140	140	20	26
103-24-080	80	80	24	30
103-24-100	100	100	24	30
103-24-120	120	120	24	30
103-24-140	140	140	24	30