



To cover tracheo- or broncho-oesophageal fistulas or in case of spontaneous or iatrogenic ruptures (such as Boerhaave syndrome) Leufen Medical offers a 80, 100 and a 120 mm stent with a diameter of 28 mm for your choice.

All oesophageal leakage stents are preloaded on the new reload delivery system. For correction of the stent position the stent can be reloaded into the delivery system during the procedure. The excellent radiopacity tantalum markers makes the stents clearly visible in the fluoroscopy.

- Stent with complete silicone covering
- Atraumatic ends
- High tensile extraction threads on both ends
- Memory effect of nitinol material – high radial force
- Radiopaque tantalum markers
- Reload delivery system

Length 700 mm, OTW

Ø 8 mm / 24 French

Accepts guide wire 0.035 inch

aixstent · oesophagus leakage · complete silicone covered

REF	Length	Length Cover	Ø Middle	Ø Ends
103-28-080	80	80	28	34
103-28-100	100	100	28	34
103-28-120	120	120	28	34