

Cleaning Brushes

Single Use

Cleaning Brushes

- Double-Sided Brushes
Efficient Cleaning of Residues
- Channel-Friendly Plastic Coated Brush Tip
Avoids Accessory Channel Scratch
- Colored Catheter
Quick Identification of Different Brush Sizes
- Stackable Dispenser Box
Hygienic and User-Friendly Storage



Cleaning Brushes

Single Use

Cleaning Brushes

Cleaning Brush – Double Sided, Individually Packed in a Dispenser Box

| Catalog No. | Brush Diameter | Catheter Length | Accessory Channel Compatibility | Unit | Product |
|---------------|-----------------------------------|-----------------|---------------------------------|------------|--|
| GAC-02-03-230 | 3 mm | 230 cm | 2.0 mm | Box of 100 |  |
| GAC-02-05-230 | 5 mm | 230 cm | 2.8 mm | Box of 100 |  |
| GAC-03-05-011 | Brush 1: 5 mm / Brush 2: 11 mm | 230 cm | 2.8 - 6.0 mm | Box of 100 |  |

Valve Cleaning Brush – Double Sided

| Catalog No. | Brush Diameter | Unit | Product |
|---------------|--------------------------------|-----------|--|
| GAC-03-05-012 | Brush 1: 5 mm / Brush 2: 12 mm | Box of 25 |  |

Cleaning Brush-Set

| Catalog No. | Cleaning Brush | | | Valve Cleaning Brush | Unit |  |
|---------------|----------------|-----------------|---------------------------------|-----------------------------------|------------|---|
| | Brush Diameter | Catheter Length | Accessory Channel Compatibility | Brush Diameter | | |
| GAC-03-05-230 | 5 mm | 230 cm | 2.8 mm | Brush 1: 5 mm / Brush 2: 12 mm | Box of 100 | |

Distributed by:

Manufactured by: