

Pediatric Gastroscopy Biteblock



Indications



Biteblocks used to facilitate the insertion of the instrument during endoscopy procedures.

**ITABLOC®**

Features



- Preloaded headstrap in laprene with differentiated elasticity.
- With tongue pad that adapts to the oral cave.
- With bilateral access for eventual aspiration
- Available without headstrap.
- Available disposable or autoclavable version.
- Packaged not sterile.
- 100 pcs. per box.

References	Diameter	Material	Colour	Headstrap
MED/106/P	16x20 mm	polypropylene	orange	•
MED/106/P/B	16x20 mm	polypropylene	orange	•
MED/106/P/A	16x20 mm	polycarbonate 	transparent	•
MED/106/P/A/B	16x20 mm	polycarbonate 	transparent	•

Specifications

- Disposable biteblock in polypropylene
- Autoclavable biteblock in polycarbonate
- Laprene headstrap with differentiated elasticity