

Pinch Clamps
Reopening Clamps



Pinch Clamps
Reopening Clamps



14065

Fits 0.059 inch -
0.118 inch OD Tubing
(1.5 mm - 3 mm)
PP

Also Available:

Part No.	Color
14064	White
14066	Red



140081

Fits 0.118 inch OD Tubing
(3 mm)
PP

Also Available:

Part No.	Material	Color
14019	ABS	White
140082	PP	White
140083	PP	Blue
140084	PP	White



140141

Fits 0.218 inch OD Tubing
(5.5 mm)
PP

Also Available:

Part No.	Color
140142	White
140143	Blue



31325

Fits 0.218 inch OD Tubing
(5.5 mm)
PP

Also Available:

Part No.	Color
31323	Red
31324	White



11089

Fits 0.125 inch -
0.15 inch OD Tubing
(3.2 mm - 3.8 mm)
PP

Also Available:

Part No.	Material	Color
14035	Acetal	White
14036	PP	White




14040

Fits 0.16 inch OD Tubing
(4 mm)
PP

Also Available:

Part No.	Color
14041	Red
14042	White



13608

Fits 0.248 inch -
0.275 inch OD Tubing
(6.3 mm - 7 mm)
POM

Also Available:

Part No.	Material	Color
13607	POM	White
13609	POM	Blue
13610	PP	Natural
13611	PP	White




14038

Fits 0.25 inch OD Tubing
(6.4 mm)
PP

Also Available:

Part No.	Color
14037	White
14039	Blue




13603

Fits 0.161 inch -
0.197 inch OD Tubing
(4.1 mm - 5 mm)
POM

Also Available:

Part No.	Material	Color
13600	POM	Natural
13601	POM	Red
13602	POM	Yellow
13604	POM	Blue
13605	PP	Natural
13606	PP	White



13020

Side Load
Fits 0.197 inch -
0.284 inch OD Tubing
(5 mm - 7.2 mm)
PP




14063

Fits 0.25 inch OD Tubing
(6.4 mm)
Zinc Plated Steel



14002

Shut Off Positive Lock
Fits 0.25 inch OD Tubing
(6.4 mm)
PP



13997

On/Off Ratchet
Fits 0.2 inch -
0.375 inch OD Tubing
(5.1 mm - 9.5 mm)
PP



11159

Lengthened Leg
Fits 0.216 inch OD Tubing
(5.4 mm)
POM

Also Available:

Part No.	Color
11160	Red
11161	White



31328

Fits 0.267 inch OD Tubing
(6.8 mm)
PP

Also Available:

Part No.	Color
31326	Red
31327	White



140073

Ratchet Style
Fits 0.268 inch OD Tubing
(6.8 mm)
PP

Also Available:

Part No.	Color
140071	Red
140072	White



13612

Fits 0.216 inch OD Tubing
(5.4 mm)
PP

Also Available:

Part No.	Color
13613	Blue




11162

Slotted
Fits 0.216 inch OD Tubing
(5.4 mm)
POM



Also Available:

Part No.	Color
11163	Blue



140322

Fits 0.268 inch OD Tubing
(6.8 mm)
PP

Also Available:


Part No.	Color
140321	White
140323	Blue



13551

Snap Locking
Fits 0.375 inch OD Tubing
(9.5 mm)
PP


Pinch Clamps Reopening Clamps



14034
Fits 0.5 inch -
0.625 inch OD Tubing
(12.7 mm - 15.9 mm)
PBT




14018
Fits 0.7 inch OD Tubing
(17.8 mm)
Acetal



11496
Fits 0.75 inch OD Tubing
(19 mm)
PBT



25659
Fits 0.75 inch OD Tubing
(19 mm)
PBT



Utilize Cable Tie (Part No. 99991) on Page 224
in order to further secure your clamp.

View a full video demonstration on qosina.com/3DPinchClamp

25659 and **99991**
Assembled with Tubing




Q9000-PC
Pinch Clamps
Sample Assortment Kit
Multiple Materials

Image Shown at 25%

QOSINA
Thousands of Stock Components


QOSINA
Thousands of Stock Components



72973
0.984 inch (25 mm)
Acetal



99940
1.2 inch (30 mm)
PP




72975
1.2 inch (30 mm)
POM



99918
1.5 inch (40 mm)
PP




99941
1.5 inch (40 mm)
PP



99942
2 inch (50 mm)
Nylon



99943
2.8 inch (70 mm)
Nylon




99915
Fits 3.5 inch (90 mm)
Nylon


Clamps, Clips
& Cable Ties

1

Pinch Clamps
Reopening Clamps



99944
4.4 inch (110 mm)
Nylon



99945
6 inch (150 mm)
Nylon



Pinch Clamps
Reopening Clamps



99914
8.625 inch (220 mm)
Sealing Length
Nylon



99988
12.5 inch (320 mm)
Sealing Length
Nylon

Pinch Clamps
Non-Reopening Clamps



13700
Securing Clamp
Fits 0.25 inch OD Tubing
(6.4 mm)
Nylon

13701
Securing Clamp
Fits 0.313 inch OD Tubing
(7.9 mm)
Nylon

13703
Securing Clamp
Fits 0.375 inch OD Tubing
(9.5 mm)
Nylon

13702
Securing Clamp
Fits 0.469 inch OD Tubing
(11.9 mm)
Nylon

72970
Nylon

72971
Nylon

72972
POM

13990
Fits 0.189 inch x
0.267 inch OD Tubing
(4.8 mm x 6.8 mm)
PP

14003
Fits 0.267 inch OD Tubing
(6.8 mm)
PP

13991
Fits 0.25 inch -
0.315 inch OD Tubing
(6.4 mm x 8 mm)
PP



51611
BioClamp™
1/2 inch - 3/4 inch
(0.5 inch - 0.75 inch;
12.7 mm - 19.1 mm)
Glass Reinforced Nylon

51612
BioClamp™
1 inch - 1 1/2 inch
(25.4 mm - 38.1 mm)
Glass Reinforced Nylon

51615
BioEnd Cap™
1/2 inch - 3/4 inch
(0.5 inch - 0.75 inch;
12.7 mm - 19.1 mm)
PP

51616
BioEnd Cap™
1 inch - 1 1/2 inch
(25.4 mm - 38.1 mm)
PP

51613
BioValve™
Fits 1/8 inch - 3/8 inch OD Tubing
(0.125 inch - 0.375 inch; 3.2 mm - 9.5 mm)
Glass Reinforced Nylon

51614
BioValve™
Fits 1/2 inch - 1 1/4 inch OD Tubing
(0.5 inch - 1.25 inch; 12.7 mm - 31.8 mm)
Glass Reinforced Nylon

Screw Clamps



14067
Fits 0.156 inch -
0.25 inch OD Tubing
(4 mm - 6.4 mm)
Acetal

Also Available:	
Part No.	Measurement
14068	Fits 0.156 inch - 0.375 inch OD Tubing (4 mm - 9.5 mm)



13576
Fits 0 inch - 0.275 inch OD Tubing
(0 mm - 6.9 mm)
PP, HDPE

Roller Clamps



14082

Fits 0.118 inch ID x
0.161 inch OD Tubing
(3 mm x 4.1 mm)
PS, HDPE



14080

Fits 0.118 inch ID x
0.161 inch OD Tubing
(3 mm x 4.1 mm)
PS, HDPE



14081

Fits 0.122 inch ID x
0.169 inch OD Tubing
(3.1 mm x 4.3 mm)
PS, HDPE, Polyolefin



33046

Fits 0.14 inch OD Tubing
(3.5 mm)
HIPS, PP




11492

Fits 0.161 inch OD Tubing
(4.1 mm)
ABS



11493
Wheel
for 11492
ABS



11494
Wheel
for 11492
ABS



X80021

Fits 0.161 inch OD Tubing
(4.1 mm)
HDPE



X80022
Wheel
for X80021
HDPE



14046

Fits 0.177 inch OD Tubing
(4.5 mm)
HDPE, ABS



14027

Fits 0.212 inch OD Tubing
(5.4 mm)
PP, HDPE



14010

Fits 0.187 inch OD Tubing
(4.7 mm)
PS, PP



14049

Fits 0.256 inch OD Tubing
(6.5 mm)
PP, ABS




QOSINA

Thousands of Stock Components


QOSINA

Thousands of Stock Components




140212

Roller Clamp Wheel
for 140211
ABS



140221

Roller Clamp Body
for 140222
Fits 0.268 inch OD Tubing
(6.8 mm)
Acetal




140222

Roller Clamp Wheel
for 140221
Acetal



140211

Roller Clamp Body
for 140212
Fits 0.268 inch OD Tubing
(6.8 mm)
ABS



140282


Roller Clamp Wheel
for 140281
ABS



33036

Roller Clamp Wheel
for 33035
ABS

Side View
with Wheel Assembled



140281

Roller Clamp Body
for 140282
Fits 0.47 inch OD Tubing
(12 mm)
ABS



33035

Roller Clamp Body
for 33036
Fits 0.47 inch OD Tubing
(12 mm)
ABS



Q9000-RC
Roller Clamps
Sample Assortment Kit
Multiple Materials



Image Shown at 25%

Roller Clamps

Clamps, Clips
& Cable Ties

Clamps, Clips
& Cable Ties

Slide Clamps



11457

Fits up to
0.06 inch OD Tubing
(0 mm - 1.5 mm)
Acetal



11099

Fits 0.07 inch -
0.12 inch OD Tubing
(1.8 mm - 3 mm)
ABS



12067

Fits 0.078 inch -
0.098 inch OD Tubing
(1.98 mm - 2.5 mm)
HDPE

Also Available:

Part No.	Color
11056	White

Also Available:

Part No.	Color
12064	White
12065	Blue
12066	Green



12069


Fits 0.078 inch -
0.098 inch OD Tubing
(1.98 mm - 2.5 mm)
HDPE

Also Available:

Part No.	Color
12068	White

QOSINA

Thousands of Stock Components




11044

Fits 0.187 inch OD Tubing
(4.7 mm)
HDPE

Also Available:

Part No.	Color
11045	Natural White



12043

Fits 0.197 inch OD Tubing
(5 mm)
PP



12053

Fits 0.197 inch OD Tubing
(5 mm)
PP



11010

Fits 0.2 inch OD Tubing
(5 mm)
HDPE



12101

Fits 0.079 inch OD Tubing
(2.01 mm)
PP



12070

Fits 0.079 inch OD Tubing
(2.01 mm)
HDPE



11026

Fits 0.079 inch -
0.161 inch OD Tubing
(2.01 mm - 4.1 mm)
PP




12071

Fits 0.098 inch
OD Tubing
(2.5 mm)
HDPE


Also Available:

Part No.	Color
12072	Red
12073	Dark Blue
12074	Light Blue



11051

Fits 0.216 inch OD Tubing
(5.4 mm)
HDPE



11038

Fits 0.25 inch OD Tubing
(6.4 mm)
HDPE

Also Available:

Part No.	Color
11039	Green
11041	Light Blue
11042	Natural White



11052

Fits 0.268 inch OD Tubing
(6.8 mm)
HDPE



13010

Fits 0.125 inch OD Tubing
(3.2 mm)
Acetal



89328

Fits 0.14 inch -
0.2 inch OD Tubing
(3.5 mm - 5 mm)
PP



12063

Fits 0.156 inch OD Tubing
(4 mm)
PP



12044

Fits 0.178 inch OD Tubing
(4.5 mm)
ABS



95867

Fits 0.275 inch OD Tubing
(7 mm)
PP



95878

Fits 0.315 inch OD Tubing
(8 mm)
ABS



95868

Fits 0.315 inch OD Tubing
(8 mm)
PP



12099

Fits 0.315 inch OD Tubing
(8 mm)
PVC

Also Available:

Part No.	Color
13014	White

Also Available:


Part No.	Color
89329	White

Slide Clamps




11029

Fits 0.315 inch OD Tubing
(8 mm)
HDPE




13015

Fits 0.375 inch OD Tubing
(9.5 mm)
ABS



11022

Fits 0.5 inch OD Tubing
(12.7 mm)
HDPE



11012

Fits 0.625 inch OD Tubing
(16 mm)
HDPE

Also Available:

Part No.	Color
11031	Green
11032	Light Blue
11033	Dark Blue
11034	Natural White

Also Available:

Part No.	Color
11021	Green
11023	Natural White
11024	Red
11025	Dark Blue


Also Available:

Part No.	Color
11013	Light Blue
11014	Natural White



Q6000-SC
Slide Clamps
Sample Assortment Kit
Multiple Materials

Image Shown at 50%




36121

0.75 inch x 21.75 inch
(19 mm x 553 mm)
10 Button Holes
Synthetic Rubber



36120

Button for
Latex-Free Strap
PS




73050

1 inch x 18 inch
(25.4 mm x 457 mm)
Polyisoprene

Also Available:

Part No.	Material	Color
73051	TPR	Blue

Flow Controllers



21292
Flow Regulator

Fits 0.118 inch ID x 0.161 inch OD Tubing
(3 mm x 4.1 mm)
ABS, NBR




17548
Flow Regulator

Fits 0.118 inch ID x 0.161 inch OD Tubing
(3 mm x 4.1 mm)
PC, Silicone



Vented
Cap Included



17549
Extension Line - Female Luer Lock
to Male Luer Lock with Flow Controller and Caps
16 inch (406 mm)
PVC, PC, Silicone, ABS, PVC, HDPE

Non-Vented
Cap Included



21293
Precision IV Flow Regulator


Fits 0.118 inch ID x 0.161 inch OD Tubing
(3 mm x 4.1 mm)
ABS, Silicone



Top View

21294
Flow Regulator

Fits 0.118 inch ID x 0.161 inch OD Tubing
(3 mm x 4.1 mm)
PC, POM, PP



Side View



54950
Oxygen Diluter
24%
PP



54951
Oxygen Diluter
28%
PP



54952
Oxygen Diluter
31%
PP




54953
Oxygen Diluter
35%
PP



54954
Oxygen Diluter
40%
PP



54955
Oxygen Diluter
50%
PP



51056
Adjustable Venturi Adapter
26% - 50%
ABS, HDPE, PP, Steel

Flow Controllers


Clamps, Clips
& Cable Ties

Cable Ties
Non-Reopening

QOSINA
Thousands of Stock Components

QOSINA
Thousands of Stock Components

Cable Ties
Reopening



11580
11 inch
(279 mm)
Nylon

99997
9 inch
(229 mm)
LDPE



99989
3 inch
(80 mm)
Nylon

99992
5.5 inch
(140 mm)
Nylon

99990
6.3 inch
(160 mm)
Nylon

99894
7 inch
(178 mm)
Nylon

99991
8 inch
(203 mm)
Nylon

99893
11 inch
(279 mm)
Nylon

99994
10 inch
(254 mm)
HDPE


99993
8 inch
(203 mm)
Nylon

999992
7 inch
(178 mm)
PE

999991
7 inch
(178 mm)
PE


99995
11 inch
(279 mm)
Nylon

Cable Ties
Reopening




- Features:**
- Insert tongue through any cell.
 - Form closed loops of various sizes with multiple wraps in either direction.
 - The unused “tail” can be cut & reused as a new tie.
 - Soft material secures items with a light, even pressure.

Retainer Bands - Latex-Free



Clips
Channel



Clips
Twist




12494
0.25 inch Diameter
(6.4 mm)
Nylon

Also Available:	
Part No.	Measurement
12495	0.37 inch Diameter (9.5 mm)
12496	0.5 inch Diameter (12.7 mm)




12497
0.24 inch Diameter
(6 mm)
PVC

Also Available:	
Part No.	Measurement
12498	0.32 inch Diameter (8 mm)
12499	0.41 inch Diameter (10.5 mm)




51268
Triple Hose Clip
0.87 inch x 3/8 inch
(22 mm x 9.5 mm)
PP



51037
Hose Clip
0.87 inch
(22 mm)
PE

NEW!




51122
Internal Gas
Sampling Line Clip
0.87 inch OD (22 mm)
PP




51496
Slide for
Cannula Tubing
0.232 inch ID (5.9 mm)
EVA

Clips
Tubing




12078
0.4 inch Diameter
(10 mm)
Nylon

Also Available:	
Part No.	Measurement
12079	0.5 inch Diameter (12.7 mm)



13570
Bed Sheet Clamp
with Catch
PP



34000
Sheet Clip, Green
HIPS

Also Available:		
Part No.	Material	Color
34001	HIPS	Blue
34002	HIPS	White
34003	HIPS	Black
34004	HIPS	Light Blue
34005	ABS	Natural


Sheet Clips



12504
0.39 inch Diameter
(10 mm)
HDPE



12500
0.46 inch Diameter
(11.6 mm)
PP



35117
Alligator Clip
Acetal



33034
Alligator Clip
with Snap Strap
Nickel Plated Steel, Vinyl

Alligator Clips

Tubing Accessories
Adhesive Backed Tubing Holders




γ 17623
1 inch Wide;
2,000 foot Roll
(25.4 mm; 610 m)
Kraft Paper, Cohesive



56498
10 FR
(0.13 inch; 3.3 mm)
PE

Peel off
and attach to nose



γ 33042
Multi-Purpose Tubing Holder
with Finger Tab
PE, Rubber



12605
Adhesive Tubing Pad/Protector
1.125 inch x 1.5 inch
(28.6 mm x 38.1 mm)
Foam, Adhesive



Tubing Accessories
Velcro® Tubing Holders



12603
Multi-Purpose Tubing Holder,
Latex-Free
Medium
HDPE/Silicone, Rayon Woven Fabric, Velcro®



12606
Tubing Positioner/Protector
1.75 inch x 1.75 inch
(44.5 mm x 44.5 mm)
Foam, Velcro®




12601
Tubing Positioner
3 inch x 3 inch
(76 mm x 76 mm)
Cotton, Velcro®




γ

11064
Tubing Hanger
0.87 inch x 1.18 inch
(22 mm x 30 mm)
PP




11050
Tubing Hanger
0.87 inch x 0.87 inch
(22 mm x 22 mm)
PP


Towel Clamp




12527
Tubing and
Wire Holder
Nylon




12526
Breathing
Circuit Clamp
Nylon




13575
Towel Clamp
ABS




56499
PP




51036
Acetal, PVC



51017
Acetal, EVA



56500
EVA, POM



12585
EVA, POM