


# RotaSnare®


Polypectomy Snares – Single Use

## RotaSnare®


- 360° Rotating Versions  
*Easier polyp snaring*
- Multifilament Loop Versions  
*Safe ectomy and coagulation*
- Monofilament Loop Versions  
*Precise and smooth cutting*
- Version with Insulated Tip  
*Minimizes risk of mucosa traumas*



### Loop Shapes



### Rotating Polypectomy Snares



#### RotaSnare® – Rotating and with Multifilament Loop

| Catalog No.   | Loop Shape              | Loop Diameter | Sheath Diameter | Working Length |
|---------------|-------------------------|---------------|-----------------|----------------|
| GPS-01-15-180 | Oval                    | 15 mm         | 2.5 mm          | 180 cm         |
| GPS-01-15-230 | Oval                    | 15 mm         | 2.5 mm          | 230 cm         |
| GPS-11-15-180 | Oval with Insulated Tip | 15 mm         | 2.5 mm          | 180 cm         |
| GPS-11-15-230 | Oval with Insulated Tip | 15 mm         | 2.5 mm          | 230 cm         |
| GPS-01-25-180 | Oval                    | 25 mm         | 2.5 mm          | 180 cm         |
| GPS-01-25-230 | Oval                    | 25 mm         | 2.5 mm          | 230 cm         |
| GPS-11-25-180 | Oval with Insulated Tip | 25 mm         | 2.5 mm          | 180 cm         |
| GPS-11-25-230 | Oval with Insulated Tip | 25 mm         | 2.5 mm          | 230 cm         |
| GPS-01-35-180 | Oval                    | 35 mm         | 2.5 mm          | 180 cm         |
| GPS-01-35-230 | Oval                    | 35 mm         | 2.5 mm          | 230 cm         |
| GPS-11-35-180 | Oval with Insulated Tip | 35 mm         | 2.5 mm          | 180 cm         |
| GPS-11-35-230 | Oval with Insulated Tip | 35 mm         | 2.5 mm          | 230 cm         |



#### RotaSnare® – Rotating and with Multifilament Loop

| Catalog No.   | Loop Shape                   | Loop Diameter | Sheath Diameter | Working Length |
|---------------|------------------------------|---------------|-----------------|----------------|
| GPS-02-25-180 | Hexagonal                    | 25 mm         | 2.5 mm          | 180 cm         |
| GPS-02-25-230 | Hexagonal                    | 25 mm         | 2.5 mm          | 230 cm         |
| GPS-12-25-180 | Hexagonal with Insulated Tip | 25 mm         | 2.5 mm          | 180 cm         |
| GPS-12-25-230 | Hexagonal with Insulated Tip | 25 mm         | 2.5 mm          | 230 cm         |

Medi-Globe polypectomy snares are directly compatible with Olympus electro-surgical generators. To verify compatibility with other electro-surgical generators please contact your Medi-Globe sales representative. Please also consult the Medi-Globe "Accessories" leaflet to see a complete list of active cords manufactured to fit most generators.

# Standard Polypectomy Snares

Oval  Oval with Insulated Tip 

## Standard with Multifilament Loop

| Catalog No.   | Loop Shape              | Loop Diameter | Sheath Diameter | Working Length |
|---------------|-------------------------|---------------|-----------------|----------------|
| GPS-21-10-180 | Oval                    | 10 mm         | 2.5 mm          | 180 cm         |
| GPS-21-10-230 | Oval                    | 10 mm         | 2.5 mm          | 230 cm         |
| GPS-21-15-180 | Oval                    | 15 mm         | 2.5 mm          | 180 cm         |
| GPS-21-15-230 | Oval                    | 15 mm         | 2.5 mm          | 230 cm         |
| GPS-31-15-180 | Oval with Insulated Tip | 15 mm         | 2.5 mm          | 180 cm         |
| GPS-31-15-230 | Oval with Insulated Tip | 15 mm         | 2.5 mm          | 230 cm         |
| GPS-21-25-180 | Oval                    | 25 mm         | 2.5 mm          | 180 cm         |
| GPS-21-25-230 | Oval                    | 25 mm         | 2.5 mm          | 230 cm         |
| GPS-21-25-300 | Oval                    | 25 mm         | 2.5 mm          | 300 cm         |
| GPS-31-25-180 | Oval with Insulated Tip | 25 mm         | 2.5 mm          | 180 cm         |
| GPS-31-25-230 | Oval with Insulated Tip | 25 mm         | 2.5 mm          | 230 cm         |
| GPS-21-35-180 | Oval                    | 35 mm         | 2.5 mm          | 180 cm         |
| GPS-21-35-230 | Oval                    | 35 mm         | 2.5 mm          | 230 cm         |
| GPS-31-35-180 | Oval with Insulated Tip | 35 mm         | 2.5 mm          | 180 cm         |
| GPS-31-35-230 | Oval with Insulated Tip | 35 mm         | 2.5 mm          | 230 cm         |
| GPS-21-50-180 | Oval                    | 50 mm         | 2.5 mm          | 180 cm         |
| GPS-21-50-230 | Oval                    | 50 mm         | 2.5 mm          | 230 cm         |

Oval Monofilament 

## Standard with Monofilament Loop

| Catalog No.   | Loop Shape | Loop Diameter | Sheath Diameter | Working Length |
|---------------|------------|---------------|-----------------|----------------|
| GPS-23-15-230 | Oval       | 15 mm         | 2.5 mm          | 230 cm         |
| GPS-23-25-230 | Oval       | 25 mm         | 2.5 mm          | 230 cm         |
| GPS-23-35-230 | Oval       | 35 mm         | 2.5 mm          | 230 cm         |

Hexagonal  Hexagonal with Insulated Tip 


## Standard with Multifilament Loop

| Catalog No.   | Loop Shape                   | Loop Diameter | Sheath Diameter | Working Length |
|---------------|------------------------------|---------------|-----------------|----------------|
| GPS-22-25-180 | Hexagonal                    | 25 mm         | 2.5 mm          | 180 cm         |
| GPS-22-25-230 | Hexagonal                    | 25 mm         | 2.5 mm          | 230 cm         |
| GPS-32-25-180 | Hexagonal with Insulated Tip | 25 mm         | 2.5 mm          | 180 cm         |
| GPS-32-25-230 | Hexagonal with Insulated Tip | 25 mm         | 2.5 mm          | 230 cm         |
| GPS-22-35-180 | Hexagonal                    | 35 mm         | 2.5 mm          | 180 cm         |
| GPS-22-35-230 | Hexagonal                    | 35 mm         | 2.5 mm          | 230 cm         |
| GPS-32-35-180 | Hexagonal with Insulated Tip | 35 mm         | 2.5 mm          | 180 cm         |
| GPS-32-35-230 | Hexagonal with Insulated Tip | 35 mm         | 2.5 mm          | 230 cm         |

Asymmetric 

## Standard with Multifilament Loop

| Catalog No.   | Loop Shape | Loop Diameter | Sheath Diameter | Working Length |
|---------------|------------|---------------|-----------------|----------------|
| GPS-41-15-230 | Asymmetric | 15 mm         | 2.5 mm          | 230 cm         |
| GPS-41-25-230 | Asymmetric | 25 mm         | 2.5 mm          | 230 cm         |



Distributed by:

Cobra Medical BV  
Stavangerweg 21-14, 9723 JC Groningen  
Postbus 5153, 9700 GD Groningen  
Nederland  
T +31 (0)50 52 80 498  
F +31 (0)50 52 80 495



[www.cobramedical.com](http://www.cobramedical.com)  
[info@cobramedical.com](mailto:info@cobramedical.com)

België / Belgique  
T 070 66 05 70  
F 070 66 05 72

Manufactured by:

Medi-Globe GmbH | Medi-Globe-Straße 1-5 | 83101 Achenmühle, Germany  
☎ +49 80 32 / 973- 379 | 📠 +49 80 32 / 973 -339 | ✉ [sales@medi-globe.de](mailto:sales@medi-globe.de)